BC 60

- (D)** Blutdruckmessgerät
Gebrauchsanleitung (2–11)
- (GB)** Blood pressure monitor
Instruction for Use (12–20)
- (F)** Tensiomètre
Mode d'emploi (21–30)
- (E)** Tensiómetro
Instrucciones para el uso (31–40)
- (I)** Misuratore di pressione
Istruzioni per l'uso (41–50)
- (TR)** Bilgisayarlı tansiyon ölçer
Kullanma Talimatı (51–59)
- (RUS)** Прибор для измерения артериального
давления на запястье
Инструкция по применению (60–70)
- (PL)** Ciśnieniomierz
Instrukcja obsługi (71–80)

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Dear Customer,

Thank you for choosing one of our products. Our name stands for high-quality, thoroughly tested products for the applications in the areas of heat, weight, blood pressure, body temperature, pulse, gentle therapy, massage and air. Please read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the information they contain.

Best regards,
Your Beurer Team

1. Getting to know your instrument

The wrist blood pressure monitor is used for non-invasive measurement and monitoring of adults' arterial blood pressure. You can use it to measure your blood pressure quickly and easily, storing the results and displaying the progression of readings together with the average.

A warning is issued for anyone suffering from cardiac arrhythmia.

The values determined are classified and graphically evaluated according to WHO guidelines.

Additionally, this blood pressure monitor includes a haemodynamic stability indicator, which will subsequently be referred to as a resting indicator. This shows whether the

circulatory system is sufficiently at rest when taking a blood pressure measurement and is therefore a more precise indicator of your resting blood pressure. Read more about this on p. 17–18.

Keep these instructions carefully for further use and also let other users have access to them.

2. Important information

Advice on use

- Always measure your blood pressure at the same time of day, in order to ensure that values are comparable.
- Relax for approx. 5–10 minutes before each measurement.
- You should wait five minutes before measuring a second time.
- The readings that you take may only be used for information purposes – they are not a substitute for a medical examination! Discuss your results with the doctor. Never use them to make medical decisions independently (e.g. regarding medication and dosage thereof).
- If circulation on one arm is restricted as a result of chronic or acute vascular disease (including vasoconstriction), the precision of any wrist measurement is impaired. In this case, use a blood pressure monitor for the upper arm.

- There may be incorrect measurements where there is cardiocirculatory disease, as well as with very low blood pressure, circulatory disorders, dysrhythmia and other preexisting diseases.
- Only use the instrument on persons whose wrist has the right circumference for the instrument.
- You may only operate the blood pressure monitor with batteries. Note that it is only possible to save data when your blood pressure monitor is receiving power. As soon as the batteries are dead, the date and time on the blood pressure monitor are lost. Any measurement results that have been stored are however retained.
- To conserve the batteries, the monitor switches off automatically if you do not press any buttons for 3 minutes.

Storage and Care

- The blood pressure monitor is made up of precision electronic components. Accuracy of readings and the instrument's service life depend on careful handling.
 - You should protect the device from impact, moisture, dirt, major temperature fluctuations and direct exposure to the sun's rays.
 - Never drop the device.
 - Do not use near strong electromagnetic fields, i.e. keep it away from any radio systems and mobile phones.

- Only ever use the cuffs provided with the monitor or original replacement cuffs. Otherwise erroneous results will be recorded.
- Do not press any buttons until the cuff is in position.
- If the instrument is not used for any length of time, we recommend removing the batteries.

Advice on batteries

- Batteries can be fatal if swallowed. You should therefore store the batteries and products where they are inaccessible to small children. If a battery has been swallowed, call a doctor immediately.
- Batteries should not be charged or reactivated with any other means, nor should they be taken apart, thrown in the fire or short-circuited.
- Remove the batteries from the instrument if they are worn out or if you are not going to use the instrument for any length of time. This prevents any damage as a result of leakage. Always replace all the batteries at the same time.
- Never use different types of battery, battery brands or batteries with different capacities. You should preferably use alkaline batteries.

Repair and disposal

- Batteries do not belong in domestic refuse. Used batteries should be disposed of at the collection points provided.

- Never open the instrument. If these instructions are not heeded, the warranty will be null and void.
- Never attempt to repair the instrument or adjust it yourself. We can no longer guarantee perfect functioning if you do.
- Repairs may only be performed by Customer Service or authorized dealers. However, always check the batteries and replace them if necessary prior to making any complaint.
- The appliance should be disposed of according to Regulation 2002/96/EC-WEEE (Waste Electrical and Electronic Equipment). In case of queries, please contact the municipal authorities responsible for waste disposal in your area.


3. Unit description

- A Display
- B Start/stop button ①
- C Memory button M1
- D Memory button M2
- E Battery compartment
- F Positioning display (heart level)

Icons in the display:


- 1 Memory (M1, M2)
- 2 WHO classification
- 3 Cardiac arrhythmia icon 
- 4 Resting indicator 
- 5 Date and time
- 6 Low battery indicator 
- 7 Pulse 
- 8 Systole/diastole 
- 9 Storage space 
- 10 Pulse Mark 


4. Getting started

Inserting battery

- Remove the battery compartment cover.
- Insert the two batteries into the battery compartment ensuring correct polarity.
- Use exclusively brand-name batteries of the type: 2x1.5 V Micro (Alkaline Type LR 03). Making absolutely sure that you insert them with the correct polarity as marked. Do not use rechargeable batteries.
- Then close the battery compartment by fitting the cover and clicking it in position.


If the low battery indicator  is flashing, then the batteries are very low. It is still possible to take a measurement, but the batteries should be replaced soon.

If the battery change  is continuously illuminated, measurement is no longer possible and you must replace all the batteries. Once batteries have been removed from the device, the time must be reset.

Used batteries do not belong in the household waste. You are legally obliged to dispose of the batteries. Dispose of them via your specialist electrical supplier or local collecting point for recyclable waste.

Note: Batteries containing pollutant substances are marked as follows: Pb = Battery contains lead, Cd = Battery contains cadmium, Hg = Battery contains mercury.

Setting date and time


- First the year flashes. Using the M2 button, you can increase the basic values to be set, using M1 you can decrease the values. Set the year and confirm the entry using the start/stop button .
- Now set the month using the M1 and M2 buttons. Confirm using the start/stop button .
- Now set the day using the M1 and M2 buttons. Confirm using the start/stop button .
- Now you can set the time, which is displayed in 24-hour mode. First set the hour, confirm using the start/stop button  and then set the minutes. After confirming with the

start/stop button,  the device switches off automatically and the time is displayed.

5. Measuring blood pressure

5.1 Positioning cuff

- Bare your left wrist, making sure that the circulation in the arm is not restricted by any clothes etc. that are too tight. Position the cuff on the inside of your wrist.
- Fasten the cuff with the Velcro fastening so that the upper edge of the monitor is positioned approx. 1 cm below the ball of your thumb.
- The cuff has to be fitted tightly around the wrist but should not constrict it.


Important: The instrument should only be operated with the original cuff.

5.2 Measuring blood pressure

- Please sit down comfortably.
- Put on the cuff as described previously and assume the position in which you want to carry out the measurement.

- Make sure to rest your arm and move it so that the device is at heart level (see image). An O.K. symbol was integrated into the device to additionally assist you. It signals the correct position of the blood pressure monitor. You are in the right position if all you can see on the display is O.K.
- Please make sure to have a 5–10 minute period of rest before measurement.
- Switch on the blood pressure monitor with the start/stop button ①.
- After the self test, during which all display elements are briefly displayed, measurement begins. The pressure is increased to 190 mmHg. If this pressure is insufficient, the device automatically pumps an extra 30 mmHG (Fuzzy Logic).
- When the measurement has been taken, the remaining air pressure dissipates very quickly. It displays pulse, systolic and diastolic blood pressure as well as the resting indicator (see Chapter 5.4).
- Now use the M1 and M2 memory buttons to select the desired user memory. If you do not select a user memory, the measurement is stored in the most recently used user memory. The relevant symbol “M1” or “M2” appears in the display.

- Using the start/stop button, ① switch off the blood pressure monitor. The measurement is then stored in the selected user memory. If you forget to turn off the device, it will switch off automatically after approx. 3 minutes. In this case too, the value is stored in the selected user memory.
- You can interrupt measurement at any time by pressing the start/stop button ①. Wait at least 5 minutes before taking another measurement!

5.3 Evaluating results

Cardiac arrhythmia:

This instrument can identify possible cardiac arrhythmia disorders during measurement and if necessary indicates the measurement with the flashing icon .

This may be an indicator for arrhythmia. Arrhythmia is a condition where the heart rhythm is abnormal as a result of defects in the bioelectrical system controlling the heart beat. The symptoms (omitted or premature heart beats, slow or excessively fast heart rate) may be caused, among other things, by heart disease, age, physical predisposition, excessive use of stimulants, stress or lack of sleep. Arrhythmia can only be ascertained through examination by your doctor. Repeat the measurement if the flashing icon  is displayed after the measurement. Please note that you should rest for 5 minutes between measurements and not talk or move during the measurement. If the icon  appears often, please contact

your doctor. Any self-diagnosis and treatment based on the test results may be dangerous.

It is vital to follow your doctor's instructions.

WHO classification:

According to WHO Guidelines/Definitions and the latest findings, the test results can be classified and evaluated according to the following chart:

Range of blood pressure values	Systolic (in mmHg)	Diastolic (in mmHg)	Measure
Grade 3: Severe hypertension	≥ 180	≥ 110	Seek medical advice
Grade 2: Moderate hypertension	160-179	100-109	Seek medical advice
Grade 1: Mild hypertension	140-159	90-99	Have it checked regularly by doctor
High-normal	130-139	85-89	Have it checked regularly by doctor
Normal	120-129	80-84	Check it yourself
Optimal	< 120	< 80	Check it yourself

Source: WHO, 1999

The bar graph in the display and the scale on the unit indicate the range of the blood pressure which has been recorded. If the values for systolic and diastolic pressure are in two different WHO ranges (e.g. systolic in the high-normal range and diastolic pressure in the normal range) the graphic WHO classification on the unit indicates the higher range (high-normal in the example described).

5.4 Resting indicator measurement (using HSD diagnostics)

The most frequent error made when measuring blood pressure is taking the measurement when not at rest (haemodynamic stability), which means that both the systolic and the diastolic blood pressures are incorrect in this case. During blood pressure measurement, the device automatically determines whether the circulatory system is sufficiently at rest or not. If there is no indication that the circulatory system is not sufficiently at rest,  (haemodynamic stability) is displayed and the measurement can be recorded as a reliable resting blood pressure value.

: **haemodynamic stability**

Measurement of the systolic and diastolic pressure is increased when the circulatory system is sufficiently at rest and is a very reliable indicator of resting blood pressure. However, if the circulatory system is not sufficiently at rest (haemodynamic instability), the symbol  is displayed. In this case, the measurement should be repeated after a period of physical and mental rest. The blood pressure measurement must be taken when the patient is physically and mentally rested, as it will be the basis for a diagnosis and regulation of the patient's medical treatment.

: **lack of haemodynamic stability**

It is very probable that the measurement of systolic and diastolic blood pressure has not been carried out at rest and has therefore distorted the measurement. Repeat the measurement after a rest and relaxation period of at least 5 minutes. Go to a sufficiently calm and comfortable spot and remain there calmly; close your eyes, breathe deeply and evenly and try to relax. If the following measurement shows a lack of stability again, repeat the measurement after another period of rest. In case of further unstable measurements, identify these measurements as they were taken when the circulatory system had not been sufficiently rested. In this case, inner unrest of the nervous system could be the cause, which cannot be cured by brief periods of rest. Devel-

oping cardiac arrhythmia can also prevent stable blood pressure measurements.

A lack of resting blood pressure can have various causes, such as physical or mental strain or distraction and speaking or experiencing heart rhythmic disturbance during the measurement.

In an overwhelming number of cases, the HSD diagnosis will give a very good guide as to whether the circulatory system is rested when taking the measurement. Certain patients suffering from cardiac arrhythmia or chronic mental conditions can remain haemodynamically unstable in the long-term, something which persists even after repeated periods of rest. The precision of the results of the resting blood pressure is reduced in these users. Like any medical measurement method, the precision of the HSD diagnosis is limited and can lead to incorrect results in some cases. The blood pressure measurements taken when the circulatory system was at rest represent particularly reliable results.

6. Calling up stored values

You can call up values from the 2 user memories and display the averages.

- Press the M1 or M2 memory buttons. Firstly, the average of all values stored in this user memory will be displayed. This is denoted by the letter "A" in the display.

- If you press the memory buttons again, further individual measurements appear with WHO, cardiac arrhythmia and HSD. These individual values are numbered from 1 to 60. The time and date are alternately displayed after 2.5 seconds.
- If there are 60 individual values in the memory at the same time, the oldest value is deleted in order to make room for the new value.
- Use the start/stop button  to turn off the device, or it will turn off automatically after approx. 30 seconds.

Tip: You can also use the 2 user memories to separately store measurements – e.g. mornings and evenings.

Deleting stored values

To delete an entire user memory, you have to select the relevant user memory. The average is displayed and denoted with the letter “A”. Press and hold memory button M1 or M2 until the display flashes and keep pressing the button until the display completely switches off.

7. Cleaning and care

- Clean your blood pressure computer carefully using only a slightly moistened cloth.
- Do not use any detergents or solvents.
- Never hold the instrument under water as otherwise liquid can penetrate and damage the instrument.
- Never place any heavy objects on the instrument.

8. Rectifying faults

Error messages can occur when

- you move or talk during the measurement (ERR 1),
- inflation takes longer than 20 seconds (ERR 2),
- device has not been inflated sufficiently despite reinflation (ERR 2),
- the inflation pressure is higher than 300 mmHg (ERR 300),
- the batteries are almost empty .

In such cases, repeat the measurement. Ensure you do not move or speak.

If necessary, re-insert or replace the batteries.

9. Specifications

Model No.	BC 60
Method of Measurement	Oscillating, non-invasive blood pressure measurement on the wrist
Range of Measurement	Pressure 0–300 mmHg, pulse 40–160 beats/minute
Accuracy of the pressure reading	systolic ± 3 mmHg / diastolic ± 3 mmHg / pulse $\pm 5\%$ of the displayed value
Measuring uncertainty	max. permissible standard deviation according to clinical testing: systolic: 8 mmHg / diastolic: 8 mmHg
Sets of Memory	2 x 60 memory places
Unit Dimension	L 72 x W 80 x H 28.7 mm
Unit Weight	Approx. 140 g

Cuff Size	For a wrist circumference of approx. 125 to 215 mm
Operating Temperature	+10 °C to +40 °C, ≤85 % RH
Storage-Temperature	-5 °C to +50 °C, ≤85 % RH
Power Supply	2 x 1,5 V Micro (Alkaline Type LR 03)
Battery Life	For 250 measurements
Accessories	Storage box, instructions for use, 2 x AAA batteries
Classification	Application part Type BF
Key to symbols:	Application part Type BF

Important! Read the operating instruction!

These specification are subject to change without notice for purpose of improvement.

- This unit is in line with European Standard EN60601-1-2 and is subject to particular precautions with regard to electromagnetic compatibility (EMC). Please note that portable and mobile HF communication systems may interfere with this instrument. For more details, please contact customer service at the address indicated.
- The instrument is in line with the EU Medical Products Directive 93/42/EC, the German medical products act and European standards EN1060-1 (non-invasive blood pressure monitors, Part 1: general requirements) and EN1060-3 (non-invasive blood pressure monitors, Part 3: Supplementary Requirements for Electromechanical Blood Pressure Measurement Systems).
- If you are using the instrument for commercial or economic purposes, you must check measurement precision regularly in accordance with the 'Operators' Ordinance on Medical Products'. Even in the case of private use, we recommend checking measurement precision every two years at the manufacturers.

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