- Ⓧ **Gebrauchsanleitung**
Blutdruckmessgerät(2-10)
- Ⓜ **Instruction for Use**
Blood pressure monitor (10-17)
- Ⓧ **Mode d'emploi**
Tensiomètre..... (18-25)
- Ⓧ **Instrucciones para el uso**
Tensiómetro (26-34)
- Ⓧ **Instruzioni per l'uso**
Misuratore di pressione (34-42)
- Ⓧ **Kullanma Talimatı**
Tansiyon Ölçme Aleti..... (43-51)
- Ⓧ **Инструкция по применению**
Прибор для измерения
артериального давления (52-60)
- Ⓧ **Instrukcja obsługi**
Ciśnieniomierz..... (60-68)

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ENGLISH

Dear Customer,

Thank you for choosing one of our products. Our name stands for high-quality, thoroughly tested products for the applications in the areas of heat, weight, blood pressure, body temperature, pulse, gentle therapy, massage and air.

Please read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the information they contain.

Best regards,
Your Beurer Team

1. Getting to know your instrument

The wrist blood pressure monitor is used for non-invasive measurement and monitoring of adults' arterial blood pressure.

You can use it to measure your blood pressure quickly and easily, storing the results and displaying the progression of readings. A warning is issued for anyone suffering from cardiac arrhythmia.

The values determined are classified and graphically evaluated according to WHO guidelines.

Keep these instructions carefully for further use and also let other users have access to them.

2. Important information

Advice on use

- Always measure your blood pressure at the same time of day, in order to ensure that values are comparable.
- Relax for approx. 5 minutes before each measurement.
- You should wait five minutes before measuring a second time.
- The readings that you take may only be used for information purposes – they are not a substitute for a medical examination! Discuss your results with the doctor. Never use them to make medical decisions independently (e.g. regarding medication and dosage thereof).
- If circulation on one arm is restricted as a result of chronic or acute vascular disease (including vasoconstriction), the precision of any wrist measurement is impaired. In this case, use a blood pressure monitor for the upper arm.
- There may be incorrect measurements where there is cardiocirculatory disease, as well as with very low blood pressure, circulatory disorders, dysrhythmia and other preexisting diseases.
- Only use the instrument on persons whose wrist has the right circumference for the instrument.
- You may only operate the blood pressure monitor with batteries. Note that it is only possible to save data when your blood pressure monitor is receiving power. As soon as the batteries are dead, the date and time on the blood pressure monitor are lost. Any measurement results that have been stored are however retained.

- The automatic switch-off function switches off the blood pressure monitor in order to preserve the batteries, if no button is pressed within one minute.

Storage and Care

- The blood pressure monitor is made up of precision electronic components. Accuracy of readings and the instrument's service life depend on careful handling.
 - You should protect the device from impact, moisture, dirt, major temperature fluctuations and direct exposure to the sun's rays.
 - Never drop the device.
 - Do not use near strong electromagnetic fields, i.e. keep it away from any radio systems and mobile phones.
 - Only ever use the cuffs provided with the monitor or original replacement cuffs. Otherwise erroneous results will be recorded.
- Do not press any buttons until the cuff is in position.
- If the instrument is not used for any length of time, we recommend removing the batteries.

Advice on batteries

- Batteries can be fatal if swallowed. You should therefore store the batteries and products where they are inaccessible to small children. If a battery has been swallowed, call a doctor immediately.
- Batteries should not be charged or reactivated with any other means, nor should they be taken apart, thrown in the fire or short-circuited.

- Remove the batteries from the instrument if they are worn out or if you are not going to use the instrument for any length of time. This prevents any damage as a result of leakage. Always replace all the batteries at the same time.
- Never use different types of battery, battery brands or batteries with different capacities. You should preferably use alkaline batteries.

i Repair and disposal

- Batteries do not belong in domestic refuse. Used batteries should be disposed of at the collection points provided.
- Never open the instrument. If these instructions are not heeded, the warranty will be null and void.
- Never attempt to repair the instrument or adjust it yourself. We can no longer guarantee perfect functioning if you do.
- Repairs may only be performed by Beurer Customer Service or authorized dealers. However, always check the batteries and replace them if necessary prior to making any complaint.
- The appliance should be disposed of according to Regulation 2002/96/EC-WEEE (Waste Electrical and Electronic Equipment). In case of queries, please contact the municipal authorities responsible for waste disposal in your area.

3. Unit description

1. Display
2. Memory button M
3. START/STOP button
4. Battery cover
5. Wrist cuff

Top of unit:

1. MODE button
2. SET button


MODE

SET

Icons in the display:

1. 'Change battery' icon
2. Error icon **ERROR**
3. Systolic pressure
4. Cardiac arrhythmia icon
5. mmHg unit

6. Icon for users 1, 2, 3
7. Diastolic pressure
8. Time and date
9. Number of memory space
10. Pulse icon ♥
11. Measured pulse
12. WHO classification

4. Prepare measurement

Inserting battery

- Remove the battery cover (left side of the unit).
- Insert two 1.5 V micro batteries (alkaline, type LR 03). Make absolutely sure that you insert the batteries with the correct polarity as marked. Do not use rechargeable batteries.
- Replace the battery cover carefully.

If the 'change battery' icon  is displayed, measurement is no longer possible and you must replace all the batteries.

Note: You will find these markings on batteries containing harmful substances:
 Pb = battery containing lead, Cd = battery containing cadmium, Hg = battery containing mercury.

Setting date and time

It is vital to set date and time. Only in this way is it possible to save and subsequently retrieve your measured values with the right date and time.

Time is shown in 24 hour format.


Note: If you hold down 'SET', you can set the values faster.

- Press 'MODE' The month starts flashing. Using the 'SET' button, set the month 1 to 12 and confirm with 'MODE'
- Set the day, hour and minute and confirm in each case with 'MODE'.

5. Measuring blood pressure

Positioning cuff

- Bare your left wrist, making sure that the circulation in the arm is not restricted by any clothes etc. that are too tight. Position the cuff on the inside of your wrist.
- Fasten the cuff with the Velcro fastening so that the upper edge of the monitor is positioned approx. 1 cm below the ball of your thumb.
- The cuff has to be fitted tightly around the wrist but should not constrict it.

Important: The instrument should only be operated with the original cuff.

Correct posture

- Rest for approx. 5 minutes before each measurement. Otherwise there may be divergences.
- You can perform the measurement either sitting or lying down. It is essential to support your arm and have it at an angle. Always make sure that the cuff is level with your heart. Otherwise there

may be serious divergences. Relax your arm and the palms of your hands.

- In order not to distort the result, it is important to keep still during the measurement and not talk.

Measuring blood pressure

- Switch on the blood pressure monitor with the 'START/STOP' button.
- Using the 'SET' button, select the user memory 1, 2 or 3. On successful measurement, results are automatically stored in the selected user memory. You can also use the different user memories to save measurement results separately, e.g. in the morning, afternoon and evening.
- Start the measurement process by pressing 'START/STOP' button. The cuff is pumped up to 195 mmHg. Cuff air pressure is released slowly. If a tendency towards high blood pressure is already detectable, the cuff is pumped up again and cuff pressure increased further. As soon as a heart rate is detected, the heart rate symbol ♥ is displayed.
- Heart rate, systolic and diastolic blood pressure are displayed.
- You can interrupt measurement at any time by pressing the 'START/STOP' button.
- The 'ERROR' icon appears if it has not been possible to perform the measurement properly. Observe the section in these instructions on error messages/troubleshooting and repeat the measurement.

Evaluating results

Cardiac arrhythmia:

This instrument can identify possible cardiac arrhythmia disorders during measurement and if necessary indicates the measurement with the flashing icon . This may be an indicator for arrhythmia. Arrhythmia is a condition where the heart rhythm is abnormal as a result of defects in the bioelectrical system controlling the heart beat. The symptoms (omitted or premature heart beats, slow or excessively fast heart rate) may

be caused, among other things, by heart disease, age, physical predisposition, excessive use of stimulants, stress or lack of sleep. Arrhythmia can only be ascertained through examination by your doctor.

Repeat the measurement if the flashing icon  is displayed after the measurement. Please note that you should rest for 5 minutes between measurements and not talk or move during the measurement. If the icon  appears often, please contact your doctor. Any self-diagnosis and treatment based on the test results may be dangerous. It is vital to follow your doctor's instructions.

WHO classification:

According to WHO Guidelines/Definitions and the latest findings, the test results can be classified and evaluated according to the following chart:

Range of blood pressure values	Systolic (in mmHg)	Diastolic (in mmHg)	Measure
Grade 3: Severe hypertension	≥ 180	≥ 110	Seek medical advice
Grade 2: Moderate hypertension	160-179	100-109	Seek medical advice
Grade 1: Mild hypertension	140-159	90-99	Have it checked regularly by doctor
High-normal	130-139	85-89	Have it checked regularly by doctor
Normal	120-129	80-84	Check it yourself
Optimal	< 120	< 80	Check it yourself

Source: WHO, 1999

The bar graph in the display and the scale on the unit indicate the range of the blood pressure which has been recorded.

If the values for systolic and diastolic pressure are in two different WHO ranges (e.g. systolic in the high-normal range and diastolic pressure in the normal range) the graphic WHO classification on the unit indicates the higher range (high-normal in the example described).

6. Saving, retrieving and deleting results

- The results of each successful measurement are stored together with date and time. With more than 30 items of measured data, the earliest items of data measured are lost.
- In order to retrieve measurement results again, press the memory button ,M'. Use the ,SET' button to select the required user memory and press the memory button ,M'. Initially, the average value of all stored results in this user memory is first displayed (display ,A'). If you continue pressing the memory button ,M', the last results in each case are measured are displayed with date, time and graphic WHO classification. With the individual stored results, warning limits and individual stored results are also displayed alternately, without pressing any other buttons.
- After pressing the memory button M and selecting the user, you can clear the memory by holding down the MODE and SET buttons simultaneously for 5 seconds.

7. Error messages/troubleshooting

In case of faults, the ,(ERROR)' message appears in the display.

9. Technical details

Model no.

Measurement method:

Measurement range

Precision of pressure reading

BC 09

Oscillating, non-invasive blood pressure measurement on the wrist

Pressure 0-300 mmHg, heart rate 40-199 beats/minute

Systolic \pm 3 mmHg / diastolic \pm 3 mmHg
pulse \pm 5% of displayed value

Error messages can occur when

1. inflation pressure is higher than 300 mmHg
2. blood pressure is unusually high or low
3. you move or talk during the measurement (in addition to ,(ERROR)', the cardiac arrhythmia icon  appears in the display)
4. the cuff tube is not inserted properly
5. the cuff is inflated for longer than 16 seconds

In the above cases, you must repeat the measurement. Make sure that the cuff tube is properly inserted and that you do not move or talk. Re-insert the batteries if necessary, or else replace them.

8. Cleaning and Care

- Clean your blood pressure computer carefully using only a slightly moistened cloth.
- Do not use any detergents or solvents.
- Never hold the instrument under water as otherwise liquid can penetrate and damage the instrument.
- Never place any heavy objects on the instrument.

Measurement uncertainty	Max. admissible standard deviation according to clinical examination.
	systolic 8 mmHg / diastolic 8 mmHg
Memory	3 x 30 storable data lines
Dimensions	L84 x W74 x H35 mm
Weight	Approx. 140 g
Cuff size	For a wrist circumference of approx. 135 to 195 mm
Admissible operating temperature	+10 °C to +40 °C, <85% relative humidity
Admissible storage temperature	-20 °C to +70 °C, <85% relative humidity
Power supply	2 x 1.5 V Micro (Alkaline Type LR 03)
Battery life	For 250 measurements
Accessories	Storage box, instructions for use, 2 x AAA batteries
Classification	Application part type BF 
Explanation of symbols	Application part type BF  Caution! Read the Instructions for Use.

- This unit is in line with European Standard EN60601-1-2 and is subject to particular precautions with regard to electromagnetic compatibility (EMC). Please note that portable and mobile HF communication systems may interfere with this instrument. For more details, please contact customer service at the address indicated.
- The instrument is in line with the EU Medical Products Directive 93/42/EC, the German medical products act and European standards EN1060-1 (non-invasive blood pressure monitors, Part 1: general requirements) and EN1060-3 (non-invasive blood pressure monitors, Part 3: Supplementary Requirements for Electromechanical Blood Pressure Measurement Systems).
- If you are using the instrument for commercial or economic purposes, you must check measurement precision regularly in accordance with the 'Operators' Ordinance on Medical Products'. Even in the case of private use, we recommend checking measurement precision every two years at the manufacturers.

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