
EIN/AUS
ON/OFF
MARCHE/ARRET

CE 0344

Typ/Type/Tipo: Maximumthermometer type FT 22/Maximum thermometer type FT 22 / Thermomètre maximum FT 22/Termómetro de máxima Tipo FT 22/Termometro di massima FT 22
Meßbereich/Range of measurement/Plage de mesure/Alcance de medición/Ambito di misurazione: 32,0 °C - 43,9 °C
Aufbewahrungstemperatur/Storage temperature/Température de stockage/Temperatura de almacenaje/Temperatura di conservazione: -10 °C - +60 °C
Normenverweis/Reference of standardisation/Norme/Referencia a normas/Normativa di riferimento: EN 12470-3; 2000

DEUTSCH

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ENGLISH

Important safety information – retain for later use.

In order to make optimum use of all the merits of the Beurer clinical thermometer, you should carefully read the operating instructions before using it, retain the instructions for later use and

also keep them accessible to other users.

- The comforter thermometer is intended only for oral measurement of body temperature in children up to 3 years of age. It is not suitable for constant sucking.
- Children should not use the thermometer without supervision.
- Before every use, check the comforter thermometer and especially the silicon nipple for signs of damage or wear.
- Damaged or worn comforter thermometers should not be used.
- Always maintain the minimum time required for measurement which is indicated by the acoustic signal. Under specific conditions, the displayed temperature may be considerably different from the actual body temperature of the child (-0.3 °C to -1.0 °C).
- Therefore, please ensure that the child does not eat or drink anything shortly before the measurement, and that the child is not crying or breathing through the mouth when the temperature is being measured.
- In case of doubt and signs of fever (temperature higher than 38.5 °C, flushed appearance, heat, perspiration, restlessness) measure the body temperature afterwards by using a conventional clinical thermometer in the anus (rectal) and discuss the measured temperatures with the physician.
- The thermometer contains sensitive electronic parts.
- Therefore, protect it from blows, bending, temperatures higher than 60 °C or direct exposure to the sun.
- Using the thermometer in strong electromagnetic fields such as next to a mobile phone can cause malfunction.
- When the thermometer is switched on, it carries out a self-test.
- Therefore, during the limited duration of use (approximately 3 years) the accuracy of measurement need not be tested.

Use:

Remove the protective covering and briefly press the ON/OFF key. A short acoustic signal will announce that the thermometer is switched on. First the thermometer will carry out a self-test for about 2 seconds. During the self-test all segments of the display will appear briefly. In case of default or a defective thermometer the default symbol "ERR" will appear (refer to "Important information").

After the self-test the measurement will start and the measuring symbol will flash "L °C". Place the comforter in the baby's mouth. Make sure the baby keeps the comforter in the mouth during the measurement and breaths calmly through the nose. The current temperature will be shown directly, e.g. 32.1 °C. The measurement takes approximately 5 minutes, which is indicated by an acoustic signal. By prolonging the length of measurement until after the acoustic signal, more accurate measurements can be obtained. If temperatures higher than 43.9 °C are measured, the symbol for too high temperature, "H" will appear. In order to increase the life of the battery, please switch off the thermometer by briefly pressing the on/off key after the measurement. Otherwise the thermometer will automatically switch off after about 10 minutes.

Memory:

The last measured value will be stored automatically. In order to display the stored temperature, press the on/off key for longer than 3 seconds. Next to the measured value the memory symbol (M) will appear. When you release the ON/OFF key the measured value will be erased and the thermometer will be ready for the next measurement.

Changing the battery/ Disposal:

The clinical thermometer contains a long-lasting silver oxide battery (IEC type SR 41) for approximately 2000 measurements. For safety reasons, the battery cannot be exchanged. When the battery symbol (▼) appears on the display, please dispose the clinical thermometer via your local dealer for electric goods or your local depot for special waste. Please dispose of the device in accordance with the directive 2002/96/EG – WEEE (Waste Electrical and Electronic Equipment). If you have any queries, please refer to the local authorities responsible for waste disposal.

Cleaning /Disinfection:

The thermometer is 100% water-proof. For the purpose of cleaning or disinfection, it can be dipped into water or disinfectant without any difficulties. Under no circumstances should the thermometer be boiled! Use one of the following disinfectants in accordance with the manufacturer's instructions. Other cleaning agents or methods could impair the function of the thermometer or damage it.

Name of the agent / Manufacturer (supplier):

70% Isopropanol, Pharmacy; Gigasept FF; Lysoformin; Aseptil.

Storage:

When you are not using the thermometer, please place the protective cap on the nipple.

Technical data:

Accuracy of measurement: ± 0.1 °C in a water bath with a temperature of 34 °C to 42 °C at an ambient temperature of 18 °C to 28 °C.

Self-test: Automatic, internal test to 37.0 °C: In case of deviations > 0.1 °C a default will be signalled by: "ERR".

Explanation of symbols: Attention – read the operating instructions.

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