## FT 55

- ⓓ **3 in 1 Ohrthermometer**  
Gebrauchsanleitung
- ⓖⓑ **3 in 1 Ear Thermometer**  
Operating instructions
- ⓕ **3 en 1 Thermomètre auriculaire**  
Mode d'emploi
- ⓔ **3 en 1 Termómetro para el oído**  
Instrucciones para el uso
- ⓐ **Termometro per orecchio 3 in 1**  
Istruzioni per l'uso
- ⓉⓇ **3'ü 1 Kulak termometresi**  
Kullanma Talimatı
- Ⓡⓞⓢ **Ушной термометр «3 в 1»**  
Инструкция по применению
- Ⓟ **Termometr douzny 3 w 1**  
Instrukcja obsługi

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## Dear customer,

We are happy that you have decided on a product from our range. Our name stands for high-quality and exhaustively tested products from the areas of warmth, gentle therapy, blood pressure, body temperature, weight, massage and air. Please read these instructions for use carefully and follow the directions.

Yours faithfully,  
The Beurer Team

## 1. Important notes

- This device is a sensitive electronic device. Please treat it carefully and do not expose it to any mechanical impacts.
- Do not expose the thermometer to direct sunlight.
- The device should be in the room in which the measurement is made for at least 30 minutes before use.
- The thermometer is NOT waterproof. For this reason, direct contact with water or other liquids should be avoided.
- Clean the measuring tip after each use with a soft cloth moistened with disinfectant.
- Check before each use that the lens is intact. If it is damaged, please contact the dealer or service address.
- The thermometer was constructed for practical use but cannot replace a visit to the doctor.
- This device is not intended for commercial or clinical use.
- The device is intended only for the purpose given in these directions for use.
- Repairs must be carried out only by authorised service agents. Otherwise the guarantee becomes void.
- There are people in whom different measurement results are obtained in the left and right ear. In order to record temperature changes, always measure in the same ear in the same person.
- If you have been lying on one ear for some time, the temperature is slightly raised. Wait a little while or measure in the other ear.
- As ear wax can affect the measurement, you should clean the ear before measuring if necessary.
- The ear thermometer may be used by children only under adult supervision. Measurement is usually possible over the age of 6 months.

In infants under 6 months, the ear passage is still very narrow so that the temperature of the eardrum can often not be recorded and the measurement result displayed is too low.

- The measurement must not be taken in an ear affected by inflammatory diseases (e.g. discharging pus or secretion), after possible ear injuries (e.g. eardrum damage) or in the healing period after operative procedures. In all of these cases, please talk to your doctor.
- This device complies with EU Directive 93/42/EC.
- This device complies with standard EN 12470-5 Medical Thermometers, Requirements for Infrared Ear Thermometers (with maximum device).
- This unit is in line with European Standard EN60601-1-2 and is subject to particular precautions with regard to electromagnetic compatibility (EMC). Please note that portable and mobile HF communication systems may interfere with this unit. For more details, please contact customer service at the address indicated.
- If you have further questions about using our devices, please contact your dealer or Customer Service.
- Please read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the information they contain.

## 2. What you should know about using this thermometer

Before using a forehead thermometer, an ear thermometer or a conventional rod thermometer, the following should be noted: The different thermometers listed above are suitable for measuring body temperature in different parts of the body: forehead thermometer: measurement only on the forehead, ear thermometer: measurement only in the ear, rod thermometer: rectal (in the back passage), axillary (under the arm) or oral measurement (in the mouth).

The temperature varies depending on the part of the body where the measurement is taken. The difference in healthy persons between different parts of the body can be between 0.2–1 °C. The approximate normal temperature ranges are as follows:

- forehead temperature – measured with a forehead thermometer: 35.8 °C to 37.6 °C.
- Ear temperature, measured with an ear thermometer: 36.0 °C to 37.8 °C.
- Rectal temperature, measured with a conventional thermometer: 36.3 °C to 37.8 °C.
- orally measured temperature – measured with a conventional thermometer: 36.0 °C to 37.4 °C.

### **Beurer-Tip:**

Temperatures measured with different thermometers cannot be compared with one another.

You should therefore tell your doctor or bear in mind if you are diagnosing yourself what thermometer you used to take your temperature and in what part of the body.

The temperature of a healthy person is also influenced by the following factors:


- The person's individual metabolism.
- Age (body temperature is higher in babies and small children and falls with increasing age. Greater temperature fluctuations occur faster and more often in children).
- Clothing.
- The outside temperature.
- The time of day (body temperature is lower in the morning and increases during the day towards evening).
- Preceding physical and, to a lesser extent, mental activity.

### Beurer-Tip:


Taking the temperature gives a measurement that provides information about a person's current body temperature. If you are uncertain about interpreting the results or if the result is abnormal (e.g. fever), you should consult your doctor. This also applies in the case of slight temperature changes if there are other symptoms of illness such as agitation, severe sweating, flushed skin, fast pulse rate, tendency to collapse etc.

## 3. Description of device

### Front side

### Rear side

## 4. Functions

This infrared thermometer measures the

- Temperature in the human ear,
- surface temperature of objects and liquids,
- ambient temperatures.

The thermometer also offers additional functionalities:

- 9 memory locations for simple plotting of changes in temperature,
- date and time, also with all stored temperatures measured,
- Optical and acoustic fever alarm at temperatures over 37.5 °C,
- switchable to °C and °F.

## 5. Operation

This thermometer is designed for use in the ear (ear passage / eardrum) usually for persons aged over 6 months. In children under 6 months, the ear passage is still very narrow so that the temperature of the eardrum can often not be recorded and the measurement result displayed is too low.

Before taking the temperature, remove the protective cap and ensure that the sensor tip and also the ear passage are clean. If there are signs of acute inflammation (discharge of pus or secretion, pain), in the case of injuries or directly after operations on the ear, the non-affected side should be chosen. Otherwise, the measurement result can be incorrect.

Use of the thermometer in different persons can be inappropriate in certain acute infectious diseases because of the possible spread of germs despite cleaning and disinfection. Discuss this with your doctor in the individual case. It must only be used without disposable covers.

There is a battery in it already so that it can be used immediately when operated for the first time.

Switch on the thermometer with the “The image shows a digital display with several icons at the top: a battery level indicator, a signal strength indicator, a Wi-Fi icon, and a 'Set' button icon. Below these icons, the display shows '24 PM 88:88' with 'Year' and 'Date' labels. At the bottom, it shows '188.8' with '°C' and '°F' labels and a small icon of a thermometer.

### 5.1 Setting time and date

With the thermometer switched on, hold the “

15

## 5.2 Batteries


The Type 3V CR2032 battery is already fitted and depending on usage has a life of around 3000 measurements. When the battery gets weak, the battery warning symbol appears . It is still possible to measure the temperature. Battery must be replaced. When the battery symbol flashes , the battery must be exchanged by sliding the battery compartment cover off to the rear.

Switch off the device, press the battery compartment lock with a sharp object and slide the battery compartment cover downwards at the same time.

Insert a new battery of the same type with the positive pole at the top. Close the battery compartment cover again.

Used batteries should not go into domestic refuse. You are legally obliged to dispose of the batteries correctly. Dispose of them through your electrical dealer or local recycling centre.

Note: you will find these symbols on batteries containing toxic substances: Pb = battery contains lead, Cd = battery contains cadmium, Hg = battery contains mercury. The battery in this device does not contain toxic substances.

## 6. During use

Always ensure that the sensor is clean and undamaged. Switch on the thermometer with the “” button.

### 6.1 Measuring in the ear

After a brief self-test and two short beeps, the thermometer is ready to measure the temperature in your ear. The infrared thermometer is in „ear thermometer“ mode. This is indicated by the symbol .

As the ear passage is slightly curved, you have to pull the ear slightly up and backwards before inserting the sensor tip. This is particularly important so that the sensor tip can be pointed directly at the eardrum. Insert the sensor tip carefully and press the “SCAN” button for 1 sec.

The end of the measuring time is signalled with a long beep. Release the “SCAN” button. You can now read off the temperature measured. The LED is also illuminated according to the test result: a green LED indicates that body temperature is within the normal range, but a red LED lights up with measured values over 37.5 °C, i.e. it is a fever alarm.

The thermometer signals with two short beeps and the cessation of flashing of the ear symbol “” that the unit is ready to continue measurement.

## 6.2 Saving the temperatures measured

The temperature last measured, i.e. only the last temperature of a series measured, is saved automatically as soon as the thermometer is switched off. There are 9 memory locations available for this.

You can call up the last temperature measured with the “⏪” button. The date and time are also shown. The LED does not support any requests for stored measured values. If you continue to hold down the “⏪”, you will enter a special measuring mode, in which the thermometer will not make a correct measurement. In this event, allow the thermometer to switch off automatically after about 1 minute and then switch it on again.

## 6.3 Surface temperatures

When you want to measure surface temperatures with this infrared thermometer, you must change to the “SCAN” mode. With the thermometer switched on (Ear thermometer mode) keep the “⏪” button depressed and at the same time press the “SCAN” button. This mode is indicated with the  symbol. If you keep the “SCAN” button depressed, the surface temperature measured is continuously displayed. You can apply the sensor tip directly to the surface to be measured or measure at close range (do not immerse in liquids).

Note that the displayed temperature is the measured and not the adapted ear surface temperature. It is not the same as the ear temperature.

## 6.4 Ambient temperature

After the measurement, the thermometer automatically switches to standby mode after about 1 minute, but only if the time has been set. This mode is identified by the  symbol and the ambient temperature is shown continuously.

If you want to use the thermometer to measure the ambient temperature, it should be positioned so that the measurement is not influenced by solar radiation or other influences such as cold draughts from air-conditioning systems. In addition to the ambient temperature, which is updated once every minute, the date and time are displayed alternately.

## 6.5 Changing the temperature unit

You can display the temperature in degrees Celsius (°C) or degrees Fahrenheit (°F). With the thermometer switched off or in standby mode, select the unit by holding down the “SCAN” button and at the same time pressing the “⏪” button. Keep the two buttons depressed until the temperature unit changes. All stored temperatures are displayed in the new temperature unit.

## 7. Trouble shooting

Error message	Problem	Solution
	Autotest in progress, device not yet ready for.	Wait until the ear symbol stops flashing.
	Large fluctuation in ambient temperature.	Store the device for at least 30 minutes in the room where the measurement will be taken.
	Ambient temperature below 10 °C or above 40 °C (<50 °F, >104 °F).	The ambient temperature must be between 10 °C and 40 °C (50 °F, 104 °F).
	Device no longer functioning correctly.	Take out battery for approx. 1 minute and insert again. If the error message is repeated, consult your dealer or customer service.
	(1) Ear thermometer mode: the temperature measured is higher than 42.2 °C (108 °F). (2) SCAN mode: the temperature measured is higher than 80 °C (176 °F).	Operate the thermometer only between the specified temperature ranges. If necessary, clean the sensor tip. In the event of a repeated error message, contact your dealer or customer service.
	(1) Ear thermometer mode: the temperature measured is lower than 34 °C (93.2 °F). (2) SCAN mode: the temperature measured is lower than -22 °C (-7.6 °F).	Operate the thermometer only between the specified temperature ranges. If necessary, clean the sensor tip. In the event of a repeated error message, contact your dealer or customer service.
	Autotest not successful.	Insert new battery.

## 8. Cleaning, storage, disposal

Clean the sensor tip after every use. Use a clean cloth or cotton bud that can be moistened with disinfectant, alcohol or warm water.

To clean the entire device, use a soft cloth slightly moistened with a mild soapy solution. Water must never get into the device. If water gets into the device, remove the battery immediately. Only use the device again when it is completely dry. Do not use any harsh cleaning agents.

Never immerse the device in water.

Always store the thermometer with the cap on in order to protect the sensor.

The device must not be stored or used at too high or low a temperature or humidity (see technical specifications), in sunlight, in association with an electrical current or in dusty locations. Otherwise, inaccuracies can occur.

If prolonged storage is intended, you should remove the battery.

The device must be tested for accuracy (calibrated) after 2 years. Please send the complete device to the dealer's or service address.

Please dispose of the device in accordance with the directive 2002/96/EG – WEEE (Waste Electrical and Electronic Equipment). If you have any queries, please refer to the local authorities responsible for waste disposal.

## 9. Technical data

**Note:** When the device is used outside the specifications, perfect functioning cannot be guaranteed.

We reserve the right to make technical changes to improve and develop the product.

Model	FT 55
Measurement range	Ear thermometer mode: 34 °C – 42.2 °C (93.2 °F – 108.0 °F) SCAN mode: -22 °C – 80 °C (-7.6 °F – 176 °F)
Laboratory measurement precision	Ear thermometer mode: ±0.2 °C (±0.4 °F) from 35.5 °C – 42 °C (95.9 °F – 107.6 °F), ±0.3 °C (±0.5 °F) < 35.5 °C > 42 °C (<95.9 °F – >107.6 °F) SCAN-Modus: ±0.3 °C (±0.5 °F) from 22 °C to 42.2 °C (71.6 °F to 108 °F) ±2 °C (±4 °F) with >42.2 °C and <22 °C (>108 °F <71.6 °F)
Clinical repeat precision	0.22 °C (children, 1 to 5 years), 0.21 °C (adults)
Measurement duration	1 to 2 seconds
Interval between 2 measurements	At least 5 seconds
Measurement units	°Celsius (°C) or °Fahrenheit (°F)
Ambient operating conditions	10 °C to 40 °C (50 °F – 104 °F) with a relative humidity of 95 % (non-condensing)

Storage conditions	-20 °C to 50 °C (-4 °F – 122 °F) with a relative humidity of up to 85 % (noncondensing)
Dimensions	34 x 145 x 28 mm
Weight	57 g including battery
Battery	1 lithium battery (type 3V CR-2032)
Memory	For 9 measurements
Explanation of symbols	Device classification type BF  Please read the instructions for use 

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