

beurer EasyFit 1.0 Installation



EASYFIT INSTALLATION

1. Place the CD-ROM supplied in the drive and choose the CD drive in Windows Explorer or on the desktop of your PC.
2. Double-click the "setup.exe" file to run it.
3. Choose your language and follow the installation instructions. Confirm with "Next" or "Install" as appropriate.
4. When you click "Finish", the software starts.
5. Enter your personal details under "New registration" in the EasyFit program and click "Save" to confirm.

DATA TRANSFER

1. Start EasyFit.
2. Go to the "Administration – Logs" menu and click the "Transfer data" button.
3. Now connect the heart rate monitor to your PC as described above. Stored data are now transferred from the heart rate monitor to the PC.

DEALING WITH PROBLEMS

Error message "Error number: 13":

Click "OK" in the message window and disconnect the USB cable from the PC. Start data transfer again. As soon as the time and date appear on the heart rate monitor, connect the heart rate monitor to the PC. Always remember to establish a new connection between the heart rate monitor and the PC for every data transfer.

If the driver has already been installed, the "Find new hardware" wizard appears:

Use the same USB port for data transfer as you used when installing the driver, or re-install the driver again (see "Driver installation") for the desired USB port.

If the heart rate monitor is already connected, "Please connect your heart rate monitor to the PC" appears in the data transfer window:

Terminate data transfer with "Cancel" and disconnect the heart rate monitor from the PC. Disconnect any scanners connected to your PC via a USB cable. Re-start the data transfer and connect the heart rate monitor to the same port as you used when installing the driver on your PC.

When scanners are connected, driver conflicts may interfere with the transfer of data.

The scanner previously disconnected from the PC can no longer be addressed:

Disconnect the heart rate monitor from your PC and restart your PC. The scanner and the heart rate monitor can now be alternately connected and used.

Data transfer can still not be successfully completed:

Uninstall the driver and then re-install it (see above). Information on uninstalling the driver can be found on our website www.beurer.de in the area "Service"/"Pulse".

SYSTEM REQUIREMENTS.

Operating system: Windows 2000/XP (SP2), CPU optimal at 1,2 GHz+, Main memory optimal at 256 MB+ RAM, Disk space optimal 200 MB+ free space - at least 100 MB, Graphic resolution starting from 1024 x 768, CD-ROM drive, USB interface.

DRIVER INSTALLATION

Basic information:

- Always use the same USB port on your PC for data transfer.
- Disconnect and connect the heart rate monitor from/to the PC every time you do a data transfer.
- Before a data transfer, disconnect any scanners connected to your PC via a USB cable.
- Note the information under "Dealing with problems".

1. First, close any applications that may be running.
2. Place the CD-ROM supplied in the drive and choose the CD drive in Windows Explorer or on the desktop of your PC.
3. Open "Drivers" subfolder in the "Software" folder.
4. Then connect the PC-Clip to the heart rate monitor and after that connect the USB plug to the PC. Be sure to connect the three pins of the clip correctly to the three contacts on the back of the heart rate monitor. The heart rate monitor display now shows "PC Link".

Windows 2000:

5. If the "Required files" window appears, close it with "Cancel".
6. Then double-click the file "DriverInstall.exe" (see step 3). After a few seconds, installation is complete. Click "OK" in the installation window.
7. Click „OK“ in the yellow window.

Windows XP:

5. If the "Find New Hardware Wizard" window appears, close it with "Cancel".
6. Then double-click the file "DriverInstall.exe" (see step 3).
7. Depending on the settings of your PC, the "Hardware installation" and/or "Software installation" may appear on the screen. These refer to the Windows logo test. Ignore these warnings and click "Continue installation" to end the installation.
8. Click „OK“ in the yellow window.

