



## THE ANTI-CELLULITE MASSAGER WITH CLINICALLY VALIDATED EFFECTS

A world first: the first clinically tested cellulite massager<sup>1)</sup>! For noticeably smoother, firmer and more toned curves<sup>2)</sup>, after just 12 weeks of use<sup>1)</sup>. The effects of the cellulite releaZer were validated in a study at the Institute proDERM

**88%** say that their skin is silky and smooth after use<sup>1)</sup>  
**94%** confirm their skin looks smoother<sup>1)</sup>

## COMBATS CELLULITE WHERE IT OCCURS



Improves visible unevenness



Promotes circulation in the layers of the skin




Reduces cellulite; dimples are less visible<sup>1)</sup>

1) Scientific study, 38 women, Institut proDERM Schenefeld/Hamburg, study no. 18.0090-53  
 2) Application results may vary



*"In the beginning I was very sceptical as I didn't believe that there really was something that could combat cellulite. I've now used the releaZer consistently for 6 weeks and already see huge improvements! I'm totally amazed and will continue with the regime<sup>2)</sup>"*




## ANTI-CELLULITE MASSAGER


Two devices, one effect. As no two women or their needs are alike, the cellulite releaZer is also available in two sizes.

### cellulite releaZer

- Reduces cellulite – scientifically proven<sup>1)</sup>
- For smooth, firm and toned curves
- Vibration massage for deep tissue toning
- 2 massage programmes/3 intensity levels
- Ergonomically shaped for easy use



**cellulite releaZer compact**  
 In pink.  
 Compact rod length: only 40 cm



**cellulite releaZer**  
 In mint green.  
 Rod length: 48 cm

Beurer GmbH,  
 P.O. Box 1427, 89004 Ulm, Soeflinger Strasse 218, 89077 Ulm, Germany  
 Tel. +49 731/39 89-0. Fax. +49 731/39 89-139.  
 export@beurer.de www.beurer.com www.beurer-healthguide.com



759.135\_08.19 Subject to errors and changes



**cellulite releaZer**  
**cellulite releaZer compact**






# Challenge your cellulite

## WHAT CAUSES CELLULITE?


Connective tissue is extremely flexible. Its collagen and elastic fibre supports are very flexible, so that they can give way well during pregnancy for example. In men, the fibres lying one above the other in a kind of net are more firmly connected to each other, meaning that there's less room for manoeuvre.

Fat cells are merciless – they exploit the flexibility of connective tissue in women and simply push themselves in between the gaps, resulting in unsightly dimples and bulges.

Cellulite is persistent: neither creams nor diets help in the long run, because they do not fight the real cause.



NORMAL SKIN



CELLULITE

## HOW DEEP VIBRATIONS COMBAT CELLULITE

The reason behind cellulite lies in the female connective tissues. Vibration massage promotes micro-circulation in the deep tissue and speeds up the reduction of fat molecules in the cells. This tones the connective tissue and makes dimples less visible. For noticeably smoother, firmer and more toned skin.<sup>2)</sup>

## THE FIRST CLINICALLY TESTED CELLULITE MASSAGER!

The effects of the cellulite releaZer were validated in a study at the Institute proDERM. Regular use of the cellulite releaZer allows you to visibly reduce persistent cellulite on the legs and buttocks.<sup>2)</sup>

The cellulite releaZer combats cellulite where it occurs. Effective vibration massage promotes the reduction of fat molecules in deep tissue, and the nub design reinforces the effects – thus resulting in legs which are more toned.<sup>2)</sup>



## FOR NOTICEABLY SMOOTH AND FIRM SKIN WITH TONED CURVES

Many women would like to have smooth thighs and a firm bottom without the unwanted dimples caused by cellulite. Now there's finally an effective solution to use at home. Combat cellulite where it occurs. Regular use allows you to visibly reduce persistent cellulite on the legs and buttocks.




## OPTIMISED TO THE LAST DETAIL!

The Beurer cellulite releaZer is an effective anti-cellulite regime for at-home use. For smooth, firm, and toned curves!

- 3 intensity levels and 2 massage programmes
- Waterproof (IPX 5) – can be used in the shower
- 4 treatment surfaces (smooth, nubs, blunt edge and sharp edge)
- Lithium-ion battery

## START YOUR CHALLENGE NOW!

The effects of the cellulite releaZer were validated in a study at the Institute proDERM. The Beurer cellulite releaZer is your effective anti-cellulite regime for at-home use. Use just every 2 days, a maximum of 15 minutes per leg.



1) Scientific study, 38 women, Institut proDERM Schenefeld/Hamburg, study no. 18.0090-53

2) Application results may vary