

SPECIAL FEATURES:

- Professional unit for use at home
- 5 intensity levels
- Integrated mirror



FC 100 Pureo Derma Peel

The FC 100 is a professional unit for use at home. The 5 intensity levels and 3 high-quality attachments mean that the device can be ideally adapted to suit your individual skin type.

- **2-way adjustable effect:** Revitalising exfoliation and an activating vacuum massage
- **Professional unit** for easy and safe use at home
- **5 intensity levels** to individually adjust to your skin type
- Illuminated open switch
- Indicator lamp
- **3 high-quality attachments** with sapphire coating
- Integrated mirror
- With 20 replacement filters
- Replacement filters available to purchase



SPECIAL FEATURES:

- Practical and portable
- Modern colour scheme



FC 76 – Microdermabrasion

The FC 76 is a practical microdermabrasion device which is also ideal for travelling.

- **2-way adjustable effect:** Revitalising exfoliation and an activating vacuum massage
- **2 intensity levels** to individually adjust to your skin type
- Indicator lamp
- **3 high-quality attachments** with sapphire coating
- With 20 replacement filters
- Replacement filters available to purchase



The renowned beauty-salon technology at home at last!

Microdermabrasion for professional facial exfoliation with the Beurer FC 100 and FC 76 can finally be achieved at home.

For a noticeably improved complexion due to refined pores, firmer and healthier-feeling skin – and all that in familiar surroundings at home without a trip to the beauty salon.



BEURER FACIAL CARE RANGE:



A. FC 45, FC 65, FC 95 B. FC 76, FC 100 C. FC 72 D. BS 89

A. Facial cleansing brushes

The water-resistant facial brushes not only gently and thoroughly cleanse facial skin, they also stimulate blood circulation at the same time. In addition to different functions and speed levels, the individual facial brushes have different brush attachments for tailored cleansing.

B. Microdermabrasion

The Pureo Derma Peel FC 100 and FC 76 devices enable professional facial peels using microdermabrasion technology. High-quality, sapphire-coated attachments are used to gently remove excess skin cells and the state-of-the-art vacuum technology noticeably improves skin texture, boosts blood circulation and stimulates cell renewal.

C. Facial sauna

The ionic facial sauna FC 72 offers the perfect basis for facial beauty treatment. The invigorating steam stimulates blood circulation in the skin and opens pores, allowing impurities to be removed more easily.

D. Illuminated cosmetics mirror

The illuminated cosmetics mirror BS 89 has a bright LED light, 5x magnification and 3x dimmer function. The mirror can be easily and conveniently attached to the wall using screws.

Interested? We would be happy to advise you!



Beurer GmbH, P.O. Box 1427, 89004 Ulm, Germany, Soefflinger Strasse 218, 89077 Ulm, Germany
Tel. +49 (0) 731/3989-0. Fax. +49 (0) 731/3989-139. export@beurer.de www.beurer.com www.beurer-healthguide.com

756.810_0816 Subject to errors and changes

beurer
beauty

Professional Microdermabrasion



Our Beauty
Our Beurer

Beautiful and well-groomed with the Beurer Beauty product range

- Professional exfoliation of your facial skin
- Activating vacuum massage
- For a noticeably improved complexion due to refined pores

Facial care products for smooth, supple skin

The skin is a reflection of your soul. A smooth, even complexion not only looks clear and nourished, it also gives you the self-assurance necessary for a confident appearance. At one time, these results could only be achieved in a beauty salon; now thanks to the Beurer beauty range you can recreate them easily in the comfort of your own home. We leave nothing to be desired.

Professional salon-quality method

Optimum facial care and treatment is achieved by combining the facial cleansing brush, microdermabrasion device and facial sauna.

STEP 1: Daily facial cleansing



Give your face a daily deep-pore cleanse with a facial cleansing brush. Then, proceed with a revitalising exfoliation or a steam treatment.

STEP 2: Microdermabrasion once a week



Carry out a revitalising exfoliation and an activating vacuum massage once a week. As you should allow your skin to rest after microdermabrasion, skip STEP 3 and finish the treatment with a skin care product.

STEP 3: Facial sauna



A steam treatment opens the pores, the skin becomes more supple and can absorb valuable active ingredients better.

STEP 4: Skin care products




To conclude the treatment, you can use a mild skin care cream or mask to pamper your skin even more.

Microdermabrasion for a refined complexion

Everyday, our skin is exposed to harmful elements such as the wind, cold, UV rays and also make-up. Microdermabrasion devices can improve skin texture through gentle exfoliation. The treatment also stimulates cell renewal for a finer and smoother complexion.

Exfoliation and vacuum massage

The revitalising exfoliation treatment gently removes excess skin particles (image 1), while the vacuum massage stimulates blood circulation and promotes cell renewal (image 2)



Attachments with sapphire coating

There are three high-quality attachments to suit different requirements and skin types. The attachment with the **fine sapphire coating** is ideal for sensitive skin. The attachment with the **rough sapphire coating** is perfect for the intensive treatment of problem areas. The **precision attachment** is ideal for the precise treatment of smaller areas such as the nose area.

Single-use filter

For hygiene reasons, you should change the filter before every use. These can be stored in the practical storage box (FC 100). Replacement filters are available to purchase.


Intensity levels

Depending on skin type and device, the user can select from up to 5 intensity levels (FC 100). This targets the skin texture more precisely and creates a more effective treatment.

Application areas

Ensure that the part of your face you are treating is taut when you move the hand-held device over the skin. You can do this by pulling the skin taut to the side or downwards with your free hand.

5 STEPS TO SUCCESS



Forehead: Move the hand-held device slowly and evenly outwards from the centre of the forehead. Make sure that you do not treat any part of your skin more than once.



Temples: Move the hand-held device outwards towards your hairline. Repeat this to the left and right of the nose. Ensure that the device is at least 1 cm away from your eyes.



Nose: Move the hand-held device downwards from the centre at the top of the bridge of your nose. Repeat this to the left and right of the nose.



Cheeks: Move the hand-held device outwards across your cheeks.



Around the mouth: Move the hand-held device along your upper lip from the centre to the left and right. Move the hand-held device upwards from your bottom lip to the left and right. Make sure that you do not touch your lips!



Application method for professional facial exfoliation

Before you start using the device, make sure that your face is clean, dry and free of any traces of products such as powder, cream or make-up.

- Place the attachment you wish to use on the tip of the hand-held device.
- Start the treatment at the lowest intensity level. Place the hand-held device on your face and move it over the skin.
- Ensure that the part of your face you are treating is taut when you move the hand-held device over the skin.
- The pressure you apply when using the device should be comfortable and not painful. Some slight reddening should be expected after use.
- If desired, you may continue your facial treatment (see "Professional salon-quality method")

PROFESSIONAL TECHNOLOGY FOR EASY APPLICATION

Is the application suitable for use on every skin type?

The microdermabrasion devices are suitable for all skin types. They can be used for both greasy and dry skin as well as combination skin.

Is the technology safe to use?

Both the Pureo-Derma Peel FC 100 and the FC 76 are based on safe microdermabrasion technology. This only removes excess skin cells from the very top layer of the skin. Other skin layers are not affected.

Is it painful to use?

The application should be gentle and pain-free. On very sensitive skin, users may experience slight burning in rare cases. If, however, you do experience pain when using the device, reduce the intensity.

How long does an average treatment last?

The length of a treatment can vary depending on the area being treated. A full facial treatment takes approx. 20 minutes.

How often should the treatment be repeated?

We recommend using microdermabrasion devices once a week and waiting at least 48 hours between treatments. You can also use the device twice a week on non-sensitive problem skin. A test must be carried out before use. Protect your skin from the sun after use.

What precautions need to be taken after treatment?

Before exposing your skin to the sun, use sun protection for the first three days after a treatment. For the first three days after a treatment DO NOT carry out any mechanical/chemical exfoliation treatments, electrolysis or facial hair removal and do not use hair removal products. You can wear make-up immediately after the treatment.