



Menstrual Relax +

TENS & heat pad



With transcutaneous electrical nerve stimulation (TENS) to provide relief from period pain



Flexible and comfortable – suitable for all body shapes



Heat function



Provides relief from period pain using TENS and heat

Provides relief from pain caused by endometriosis (painful chronic condition affecting women)

Larger shape and heat surface for added warmth

For use at home and on the move

2 in 1: heat and TENS function

Medication-free pain relief

Soft, flexible and discreet - Fits to your body shape

1 heat programme (43°C), 15 TENS intensity levels

Heat and TENS can be selected separately

LED display for the set function

Up to 6 applications per battery charge (6 x 20 min)

Application duration: 20 min

Automatic switch-off and safety switch-off function

Incl. 4 self-adhesive gel pads and storage bag

Reusable

Indicator shows when battery state is low

Rechargeable (incl. USB cable)

Medical device

Product dimensions (LxWxH): approx. 235 x 115 x 20 mm

Product weight: approx. 82 g

5-year warranty

Pack: 6 / Shipping carton: 8 x 6

EAN no.: 4211125100599

Item no.: 100.59

