




EM 50 Menstrual Relax TENS & heat pad

 german
engineering




With transcutaneous electrical nerve stimulation (TENS) to provide relief from period pain



15 TENS intensity levels



Heat function



Provides relief from period pain using TENS and heat

Provides relief from pain caused by endometriosis (painful chronic condition affecting women)

For use at home and on the move

2 in 1: heat and TENS function

Medication-free pain relief

Soft and flexible

Comforting and discreet

Fits to your body shape

1 heat programme (43°C), 15 TENS intensity levels

Heat and TENS can be selected separately

LED display for the set function

Up to 6 applications per battery charge (6 x 20 min)

Application duration: 20 min

Automatic switch-off and safety switch-off function

Incl. 2 self-adhesive gel pads

Reusable

Indicator shows when battery state is low

Rechargeable (incl. USB cable)

Medical device

Product dimensions (LxWxH): approx. 95 x 193 x 15 mm

Product weight: approx. 47 g

5-year warranty

Pack: 6 / Shipping carton: 12 x 6

EAN no.: 4211125648541

Item no.: 648.54