




BM 64

beurer
medical

BM 64 - With atrial fibrillation (AFib) detection for stroke prevention



Advanced arrhythmia detection:
Detects atrial fibrillation for
stroke prevention



Triple measurement is recommended
by doctors for even more accurate
measurements



Inflation technology: Ultimate measure-
ment comfort thanks to gentle build-up
of pressure and quick measuring time



Fully automatic blood pressure and pulse measurement on the upper arm

Arrhythmia detection: detects and warns of possible cardiac arrhythmias

Patented resting indicator* for exact measurement results

Cuff position check: For ensuring the cuff is correctly positioned on the upper arm

Risk Indicator: Classification of measurements with coloured scale

Two user memories for 120 measured values each

Easy to use via 4 sensor touch buttons

Monitor changes in your blood pressure values at all times with the "beurer HealthManager Pro" app

Easy-to-read black XL display with backlight

Calculation of the average value of all saved measured values and morning and evening blood pressure for the last 7 days

Automatic pressure preselection and deflation

Medical device

Date and time /automatic switch-off

Battery changing indicator

With premium cuff for upper arm circumferences from 22 - 42 cm

Innovative networking between smartphone and blood pressure monitor – keep track of your measurements at home and on the move

Incl. 4 x 1,5V AA batteries, USB cable and storage pouch

Product measurements (LxWxH): 143 mm x 105 mm x 60 mm

Product weight: approx. 450 g (without batteries, with cuff)

5 year guarantee

Sales unit: 4 / Shipping carton: 4 x 3

EAN no.: 4 211125 65400 9

Item no.: 654.00

* Patent EP 1 673 009

654.00_BM64_2023-06-16_01_DS_EN_BEU Irrtum und Änderungen vorbehalten

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Beurer GmbH is under license. Other trademarks and trade names are those of their respective owners.

Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries.

Google Play and the Google Play logo are trademarks of Google LLC.

Android is a trademark of Google LLC.