



PelviStrong

Pelvic Floor Training Aids



Made in Germany



Regular exercise can strengthen the pelvic floor



Three cones of the same size
with different weights



3 tampon-shaped cones with different weight: 28 g, 52 g, 71 g

Incl. exercises with training instructions

Comfortable, smooth surface offers easy cleaning after use

Produktmaße (LxB): 51 x 23 mm

Produktgewicht:

turquoise: 28 g

fir green: 52 g

midnight blue: 71 g

VE: 6 / Transportkarton: 4 x 6

EAN-Nr.: 421112 579072 1

Art.-Nr.: 790.72